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Test Name: NN **D-Dimer** Test Code: IL_ **DDim-d** **Measure** **Sample Reagent** **Ranges Limits** **Calibration**

Sample Type: **Serum**

Reaction Cycle: ☒ Standard ☐ Extended

Reporting Unit: ☒ МКг FEU /мл * Decimal Points: 2 ☐ User Define Conv. Factor: 0.000

Methodology: Type: ☐ End Point ☒ Rate Measuring Point: 19 - ☒ 30 Photometric: ☒ 1 Wavelength ☐ 2 Wavelength Primary: 570 Secondary: Correction Constant: Slope: 1.000 Intercept: 0.000

Test Name: NN **D-Dimer** Test Code: IL_ **DDim-d** **Measure** **Sample Reagent** **Ranges Limits** **Calibration**

Sample Volume:

	Sample	Dilution	
	Volume	Sample Vol.	Diluent Vol.
1	4.0	0.0	0.0
2	0.0	0.0	0.0
3	0.0	0.0	0.0
4	0.0	0.0	0.0

 Diluent: Diluent Warning Limit: 0 Tests

Reagent Volume: R1: **DDim-d1** R2: **DDim-d2**

	Volume	Diluent Vol.	stirring	Warning Limit(tests)	Stability(days)
R1	120	0.0	<input checked="" type="checkbox"/>	20	No Control
R2	60	0.0	<input checked="" type="checkbox"/>	20	No Control

Sampling Condition:

	Condition No.	1	2	3	4
<input type="checkbox"/>	First Run	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Samp.Vol.Reduction	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Below N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Above N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Panic L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Panic H	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	User Range L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	User Range H	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Noise	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Prozone	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	HIGH!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	ABS!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>

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Ranges		Limits																																												
<table border="1"><thead><tr><th colspan="3">Normal Range</th></tr><tr><th></th><th>Lower</th><th>Upper</th></tr></thead><tbody><tr><td>Male</td><td>0.0</td><td>0.5</td></tr><tr><td>Female</td><td>0.0</td><td>0.5</td></tr><tr><td>Other</td><td>0.0</td><td>0.5</td></tr></tbody></table>		Normal Range				Lower	Upper	Male	0.0	0.5	Female	0.0	0.5	Other	0.0	0.5	<table border="1"><thead><tr><th>Reaction Slope</th><th>Absorbance Limit</th></tr></thead><tbody><tr><td><input type="radio"/> Negative <input checked="" type="radio"/> Positive</td><td><input checked="" type="radio"/> Above <input type="radio"/> Below</td></tr><tr><td><table border="1"><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>100</td><td>%</td></tr></thead></table></td><td><table border="1"><thead><tr><td>2600.0</td><td>mAbs</td></tr></thead></table></td></tr><tr><td colspan="2"><table border="1"><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td><table border="1"><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table></td><td></td></tr><tr><td colspan="2"><table border="1"><thead><tr><th>Judge Point</th></tr><tr><td>0</td></tr></thead></table></td></tr></thead></table></td></tr></tbody></table>		Reaction Slope	Absorbance Limit	<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below	<table border="1"><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>100</td><td>%</td></tr></thead></table>	Non-Linear Limit		100	%	<table border="1"><thead><tr><td>2600.0</td><td>mAbs</td></tr></thead></table>	2600.0	mAbs	<table border="1"><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td><table border="1"><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table></td><td></td></tr><tr><td colspan="2"><table border="1"><thead><tr><th>Judge Point</th></tr><tr><td>0</td></tr></thead></table></td></tr></thead></table>		Prozone Limit		<input type="radio"/> Above <input checked="" type="radio"/> Below		<table border="1"><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table>	Limit	Equation	0.0	none		<table border="1"><thead><tr><th>Judge Point</th></tr><tr><td>0</td></tr></thead></table>		Judge Point	0
Normal Range																																														
	Lower	Upper																																												
Male	0.0	0.5																																												
Female	0.0	0.5																																												
Other	0.0	0.5																																												
Reaction Slope	Absorbance Limit																																													
<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below																																													
<table border="1"><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>100</td><td>%</td></tr></thead></table>	Non-Linear Limit		100	%	<table border="1"><thead><tr><td>2600.0</td><td>mAbs</td></tr></thead></table>	2600.0	mAbs																																							
Non-Linear Limit																																														
100	%																																													
2600.0	mAbs																																													
<table border="1"><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td><table border="1"><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table></td><td></td></tr><tr><td colspan="2"><table border="1"><thead><tr><th>Judge Point</th></tr><tr><td>0</td></tr></thead></table></td></tr></thead></table>		Prozone Limit		<input type="radio"/> Above <input checked="" type="radio"/> Below		<table border="1"><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table>	Limit	Equation	0.0	none		<table border="1"><thead><tr><th>Judge Point</th></tr><tr><td>0</td></tr></thead></table>		Judge Point	0																															
Prozone Limit																																														
<input type="radio"/> Above <input checked="" type="radio"/> Below																																														
<table border="1"><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table>	Limit	Equation	0.0	none																																										
Limit	Equation																																													
0.0	none																																													
<table border="1"><thead><tr><th>Judge Point</th></tr><tr><td>0</td></tr></thead></table>		Judge Point	0																																											
Judge Point																																														
0																																														

User Range	
Lower	Upper

ValidRange	
Lower	Upper
0.1	8.7

Qualitative
☐ On ☒ Off

Test Name		Test Code		Measure	Sample Reagent	Ranges Limits	Calibration
NN	D-Dimer	IL_	DDim-d				

Calibration	Condition	Calibration/R-Blank Limit																																																								
<table border="1"><thead><tr><th>Method</th><th>Curve Type</th></tr></thead><tbody><tr><td>M-point</td><td>5P Log</td></tr></tbody></table> <table border="1"><thead><tr><th>Repeats</th><th>Stability</th></tr></thead><tbody><tr><td>2</td><td></td></tr><tr><td>No</td><td></td></tr><tr><td>6</td><td></td></tr></tbody></table> <table border="1"><thead><tr><th colspan="3">Calibrator</th></tr><tr><th></th><th>Calibrator</th><th>Conc</th></tr></thead><tbody><tr><td>1</td><td>D-Diluent</td><td>0.0</td></tr><tr><td>2</td><td>TC-DDim-1</td><td>**</td></tr><tr><td>3</td><td>TC-DDim-2</td><td>**</td></tr><tr><td>4</td><td>TC-DDim-3</td><td>**</td></tr><tr><td>5</td><td>TC-DDim-4</td><td>**</td></tr><tr><td>6</td><td>TC-DDim-5</td><td>**</td></tr></tbody></table>	Method	Curve Type	M-point	5P Log	Repeats	Stability	2		No		6		Calibrator				Calibrator	Conc	1	D-Diluent	0.0	2	TC-DDim-1	**	3	TC-DDim-2	**	4	TC-DDim-3	**	5	TC-DDim-4	**	6	TC-DDim-5	**	<div><input type="checkbox"/> Reagent Blank</div> <div><input type="checkbox"/> Auto Reagent Blank by Bottle</div> <div><input type="checkbox"/> Ask for calibration when reagent lot changes</div>	<table border="1"><thead><tr><th colspan="2">R-Blank Limit</th></tr><tr><td>2500.0</td><td>mAbs</td></tr></thead></table> <table border="1"><thead><tr><th colspan="2">Cal Reps Range</th></tr><tr><td>20.0</td><td>%</td></tr></thead></table> <table border="1"><thead><tr><th colspan="2">Min Cal Reps</th></tr><tr><td>100.0</td><td>mAbs</td></tr></thead></table> <table border="1"><thead><tr><th colspan="2">Factor Change</th></tr><tr><td>35</td><td>%</td></tr></thead></table> <table border="1"><thead><tr><th colspan="2">M-Point Curve Fit</th></tr><tr><td>0.0</td><td>%</td></tr></thead></table>	R-Blank Limit		2500.0	mAbs	Cal Reps Range		20.0	%	Min Cal Reps		100.0	mAbs	Factor Change		35	%	M-Point Curve Fit		0.0	%
Method	Curve Type																																																									
M-point	5P Log																																																									
Repeats	Stability																																																									
2																																																										
No																																																										
6																																																										
Calibrator																																																										
	Calibrator	Conc																																																								
1	D-Diluent	0.0																																																								
2	TC-DDim-1	**																																																								
3	TC-DDim-2	**																																																								
4	TC-DDim-3	**																																																								
5	TC-DDim-4	**																																																								
6	TC-DDim-5	**																																																								
R-Blank Limit																																																										
2500.0	mAbs																																																									
Cal Reps Range																																																										
20.0	%																																																									
Min Cal Reps																																																										
100.0	mAbs																																																									
Factor Change																																																										
35	%																																																									
M-Point Curve Fit																																																										
0.0	%																																																									

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