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Test Name: NN **EtOH** Test Code: IL_ **EtOH-d** **Measure** Sample Reagent Ranges Limits Calibration

Sample Type
Serum

Reaction Cycle
☒ Standard ☐ Extended

Reporting Unit
☒ г/л* Decimal Points: **2**
☐ User Define
 Conv. Factor: **0.000**

Methodology
 Type: ☒ End Point ☐ Rate
 Measuring Point: **16** - ☒ **33**
 Photometric: ☐ 1 Wavelength ☒ 2 Wavelength
 Primary: **375** Secondary: **546**
 Correction Constant:
 Slope: **1.000** Intercept: **0.000**

Test Name: NN **EtOH** Test Code: IL_ **EtOH-d** **Measure** **Sample Reagent** Ranges Limits Calibration

Sample Volume

	Sample	Dilution	
	Volume	Sample Vol.	Diluent Vol.
1	2.0	0.0	0.0
2	2.0	0.0	0.0
3	4.0	0.0	0.0
4	0.0	0.0	0.0

Diluent

Diluent Warning Limit
0 Tests

Reagent Volume
 R1: **EtOH-d1** R2: **EtOH-d2**

	Volume	Diluent Vol.	stirring	Warning Limit(tests)	Stability(days)
R1	200	0.0	<input checked="" type="checkbox"/>	10	No Control
R2	50	0.0	<input checked="" type="checkbox"/>	10	No Control

Sampling Condition

	Condition No.	1	2	3	4
<input type="checkbox"/>	First Run	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Samp.Vol.Reduction	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Below N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Above N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Panic L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Panic H	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	User Range L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	User Range H	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Noise	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Prozone	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	HIGH!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	ABS!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>

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Ranges		Limits																																		
<table border="1"><thead><tr><th colspan="3">Normal Range</th></tr><tr><th></th><th>Lower</th><th>Upper</th></tr></thead><tbody><tr><td>Male</td><td>0.0</td><td>0.05</td></tr><tr><td>Female</td><td>0.0</td><td>0.05</td></tr><tr><td>Other</td><td>0.0</td><td>0.05</td></tr></tbody></table>		Normal Range				Lower	Upper	Male	0.0	0.05	Female	0.0	0.05	Other	0.0	0.05	<table border="1"><thead><tr><th>Reaction Slope</th><th>Absorbance Limit</th></tr></thead><tbody><tr><td><input type="radio"/> Negative <input checked="" type="radio"/> Positive</td><td><input checked="" type="radio"/> Above <input type="radio"/> Below</td></tr><tr><td><table border="1"><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>0</td><td>%</td></tr></thead><tbody></tbody></table></td><td><table border="1"><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td>Limit: 0.0</td><td>Equation: none</td></tr><tr><td colspan="2">Judge Point: 0</td></tr></thead></table></td></tr></tbody></table>		Reaction Slope	Absorbance Limit	<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below	<table border="1"><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>0</td><td>%</td></tr></thead><tbody></tbody></table>	Non-Linear Limit		0	%	<table border="1"><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td>Limit: 0.0</td><td>Equation: none</td></tr><tr><td colspan="2">Judge Point: 0</td></tr></thead></table>	Prozone Limit		<input type="radio"/> Above <input checked="" type="radio"/> Below		Limit: 0.0	Equation: none	Judge Point: 0	
Normal Range																																				
	Lower	Upper																																		
Male	0.0	0.05																																		
Female	0.0	0.05																																		
Other	0.0	0.05																																		
Reaction Slope	Absorbance Limit																																			
<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below																																			
<table border="1"><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>0</td><td>%</td></tr></thead><tbody></tbody></table>	Non-Linear Limit		0	%	<table border="1"><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td>Limit: 0.0</td><td>Equation: none</td></tr><tr><td colspan="2">Judge Point: 0</td></tr></thead></table>	Prozone Limit		<input type="radio"/> Above <input checked="" type="radio"/> Below		Limit: 0.0	Equation: none	Judge Point: 0																								
Non-Linear Limit																																				
0	%																																			
Prozone Limit																																				
<input type="radio"/> Above <input checked="" type="radio"/> Below																																				
Limit: 0.0	Equation: none																																			
Judge Point: 0																																				

	User Range			------------	-------		Lower	Upper		0	3.50				
	ValidRange			------------	-------		Lower	Upper		0	3.50				
Qualitative															
☐ On ☒ Off															

Test Name		Test Code		Measure		Sample Reagent		Ranges Limits		Calibration																																		
NN	EtOH	IL	EtOH-d																																									
Calibration				Condition				Calibration/R-Blank Limit																																				
<table border="1"><thead><tr><th colspan="2">Method</th></tr><tr><td>1-point</td><td></td></tr></thead><tbody></tbody></table>				Method		1-point		<table border="1"><thead><tr><th colspan="2">Curve Type</th></tr><tr><td>Linear</td><td></td></tr></thead><tbody></tbody></table>				Curve Type		Linear		<table border="1"><thead><tr><th colspan="2">R-Blank Limit</th></tr><tr><td>2500.0</td><td>mAbs</td></tr></thead><tbody></tbody></table>				R-Blank Limit		2500.0	mAbs																					
Method																																												
1-point																																												
Curve Type																																												
Linear																																												
R-Blank Limit																																												
2500.0	mAbs																																											
<table border="1"><thead><tr><th colspan="2">Repeats</th></tr><tr><td>2</td><td></td></tr></thead><tbody></tbody></table>				Repeats		2		<table border="1"><thead><tr><th colspan="2">Stability</th></tr><tr><td></td><td>days</td></tr></thead><tbody></tbody></table>				Stability			days	<table border="1"><thead><tr><th colspan="2">Cal Reps Range</th></tr><tr><td>20.0</td><td>%</td></tr></thead><tbody></tbody></table>				Cal Reps Range		20.0	%																					
Repeats																																												
2																																												
Stability																																												
	days																																											
Cal Reps Range																																												
20.0	%																																											
<table border="1"><thead><tr><th colspan="2">No</th></tr><tr><td></td><td></td></tr></thead><tbody></tbody></table>				No				<table border="1"><thead><tr><th colspan="2">Min Cal Reps</th></tr><tr><td>100.0</td><td>mAbs</td></tr></thead><tbody></tbody></table>				Min Cal Reps		100.0	mAbs	<table border="1"><thead><tr><th colspan="2">Factor Change</th></tr><tr><td>20</td><td>%</td></tr></thead><tbody></tbody></table>				Factor Change		20	%																					
No																																												
Min Cal Reps																																												
100.0	mAbs																																											
Factor Change																																												
20	%																																											
<table border="1"><thead><tr><th colspan="3">Calibrator</th></tr><tr><th></th><th>Calibrator</th><th>Conc</th></tr></thead><tbody><tr><td>1</td><td>Std_EtOH</td><td>**</td></tr><tr><td>2</td><td></td><td></td></tr><tr><td>3</td><td></td><td></td></tr><tr><td>4</td><td></td><td></td></tr><tr><td>5</td><td></td><td></td></tr></tbody></table>				Calibrator				Calibrator	Conc	1	Std_EtOH	**	2			3			4			5			<table border="1"><thead><tr><th colspan="2">M-Point Curve Fit</th></tr><tr><td>0.0</td><td>%</td></tr></thead><tbody></tbody></table>				M-Point Curve Fit		0.0	%	<table border="1"><thead><tr><th colspan="2">Reagent Blank</th></tr><tr><td><input checked="" type="checkbox"/></td><td></td></tr></thead><tbody><tr><td><input type="checkbox"/></td><td>Auto Reagent Blank by Bottle</td></tr><tr><td><input type="checkbox"/></td><td>Ask for calibration when reagent lot changes</td></tr></tbody></table>				Reagent Blank		<input checked="" type="checkbox"/>		<input type="checkbox"/>	Auto Reagent Blank by Bottle	<input type="checkbox"/>	Ask for calibration when reagent lot changes
Calibrator																																												
	Calibrator	Conc																																										
1	Std_EtOH	**																																										
2																																												
3																																												
4																																												
5																																												
M-Point Curve Fit																																												
0.0	%																																											
Reagent Blank																																												
<input checked="" type="checkbox"/>																																												
<input type="checkbox"/>	Auto Reagent Blank by Bottle																																											
<input type="checkbox"/>	Ask for calibration when reagent lot changes																																											

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