

## Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: CRPu Type: Photometric Version: 1  
 Number: \*\*\* Assay availability: Enabled Date: 18.02.2020  
 Run controls for onboard reagents by: Lot Time: 11:01:08  
 Operator: ADMIN

Reaction definition Reagent / Sample Validity checks

Reaction mode: Rate up

Primary Secondary Read times  
 Wavelength: 505 / None Main: 19 -- 27  
 Last required read: 33 Flex: --  
 Absorbance range: -- Color correction: --  
 Sample blank type: None

Assay of

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: CRPu Type: Photometric Version: 1  
 Number: \*\*\* Assay availability: Enabled Date: 20.02.2020  
 Run controls for onboard reagents by: Lot Time: 17:14:07  
 Operator: FSE

Reaction definition Reagent / Sample Validity checks

Reagent: CRPu  
 Diluent: Saline  
 Diluent dispense mode: Type 0

R1 R2  
 Reagent volume: 80 80  
 Water volume: --  
 Dispense mode: Type 2 Type 2

Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor	Default dilution
Standart	1.6				1:1.00	<input checked="" type="radio"/>
1:5	20.0	1.6	80.0		1:5.00	<input type="radio"/>
2x	3.2				2:1.00	<input type="radio"/>

Assay of

Done Cancel

?

## Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: CRPu Assay number: \*\*\* Date: 20.02.2020  
Calibration method: Spline Time: 17:14  
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator set: Calibrator level: Concentration:

TruCalCRPu Blank: Water 0.0000

Cal 1: TruCalCRPu1  
Cal 2: TruCalCRPu2  
Cal 3: TruCalCRPu3  
Cal 4: TruCalCRPu4  
Cal 5: TruCalCRPu5  
Cal 6: None

Replicates: 2 [Range 1 - 3]

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: CRPu Assay number: \*\*\* Date: 20.02.2020  
Calibration method: Spline Time: 17:14  
Operator: FSE

Calibrators Volumes Intervals Validity checks


Calibrator:

Calibrator level	Sample	Diluted sample	Diluent	Water
Blank: Water	1.6			
Cal 1: TruCalCRPu1	1.6			
Cal 2: TruCalCRPu2	1.6			
Cal 3: TruCalCRPu3	1.6			
Cal 4: TruCalCRPu4	1.6			
Cal 5: TruCalCRPu5	1.6			
Cal 6: None				

Done Cancel

?

## Параметры для ввода в программу анализатора Architect c8000 / c4000



Единицы измерения - мг/л; знаков после запятой 2 (**Configure result units**).

Диапазон метода: 0.2-350 мг/л (**Configure assay parameters\ Results**).

Референсный интервал: 0-5 мг/л (**Configure assay parameters\ Results**).  
При наличии ЛИС, не обязательно.

Диапазон нормальных значений указан в соответствии с рекомендациями производителя реагентов.

Калибратор TruCal CRP-U, набор из 5 уровней, нулевой - вода (или физраствор).

Контроль по TruLab CRP уровень 1 и уровень 2 или TruLab Protein уровень 1 и уровень 2.