

Имуноглобулин А, Ig-A

ImmunoTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: IgA Type: Photometric Version: 1
Number: *** Assay availability: Enabled Date: 18.02.2020
Run controls for onboard reagents by: Lot Time: 11:01:08
Operator: ADMIN

Reaction definition Reagent / Sample Validity checks

Reaction mode: End up

Primary Secondary Read times
Wavelength: 572 / None Main: 30 -- 31
Last required read: 33
Absorbance range: 0.000 -- 2.500 Color correction: --
Sample blank type: Self Blank: 14 -- 16

Assay 13 of 54

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: IgA Type: Photometric Version: 1
Number: *** Assay availability: Enabled Date: 20.02.2020
Run controls for onboard reagents by: Lot Time: 17:14:07
Operator: FSE

Reaction definition Reagent / Sample Validity checks

Reagent: IgA
Diluent: Saline
Diluent dispense mode: Type 0

R1 R2
Reagent volume: 190 38
Water volume: --
Dispense mode: Type 2 Type 2

Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor	Default dilution
Standart	1.5				1:1.00	<input checked="" type="radio"/>
1:4	20.0	1.5	80		1:4.00	<input type="radio"/>
						<input type="radio"/>


Assay of

Done Cancel

?

Имуноглобулин А, Ig-A

ImmunoTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: IgA Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator set: TruCalProt
Replicates: 2 [Range 1 - 3]

Calibrator level: Blank: Water
Cal 1: TruCalProt1
Cal 2: TruCalProt2
Cal 3: TruCalProt3
Cal 4: TruCalProt4
Cal 5: TruCalProt5
Cal 6: None

Concentration: 0.0000

Assay of

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: IgA Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator: TruCalProt

Calibrator level	Sample	Diluted sample	Diluent	Water
Blank: Water	1.5			
Cal 1: TruCalProt1	1.5			
Cal 2: TruCalProt2	1.5			
Cal 3: TruCalProt3	1.5			
Cal 4: TruCalProt4	1.5			
Cal 5: TruCalProt 5	1.5			
Cal 6: None				

Assay of

Done Cancel

?

Имуноглобулин А, Ig-A

ImmunoTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

☐ General ☒ Calibration ☐ SmartWash ☐ Results ☐ Interpretation

Assay: IgA Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

☐ Calibrators ☐ Volumes ☒ Intervals ☐ Validity checks

Calibration intervals:

Full interval: 0 (hours)
Adjust interval: 0 (hours)

Calibration type:

Adjust type: 2 Point
Adjust level: Cal 1
Default ordering type: Full

Assay of

Done
Cancel

?

Единицы измерения - г/л; знаков после запятой 2 (**Configure result units**).

Диапазон метода: 0.3 - 9.0 г/л (**Configure assay parameters\ Results**).

Рефереренсный интервал: 0.7-4.0 г/л (**Configure assay parameters\ Results**).

При наличии ЛИС, не обязательно.

Диапазон нормальных значений указан в соответствии с рекомендациями производителя реагентов.

Калибратор TruCal Protein, набор из 5 уровней, нулевой - вода (или физраствор).

Контроль по TruLab Protein уровень 1 и уровень 2.