

# Имуноглобулин М, Ig-M

ИммуноТД (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect с8000 / с4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: IgM Type: Photometric Version: 1  
Number: \*\*\* Assay availability: Enabled Date: 18.02.2020  
Run controls for onboard reagents by: Lot Time: 11:01:08  
Operator: ADMIN

Reaction definition Reagent / Sample Validity checks

Reaction mode: End up

Primary Secondary Read times  
Wavelength: 404 / 804 Main: 30 -- 31  
Last required read: 33  
Absorbance range: 0.000 -- 2.500 Color correction: --  
Sample blank type: Self Blank: 14 -- 16

Assay 13 of 54

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: IgM Type: Photometric Version: 1  
Number: \*\*\* Assay availability: Enabled Date: 20.02.2020  
Run controls for onboard reagents by: Lot Time: 17:14:07  
Operator: FSE

Reaction definition Reagent / Sample Validity checks

Reagent: IgM  
Diluent: Saline  
Diluent dispense mode: Type 0

R1 R2  
Reagent volume: 190 38  
Water volume: --  
Dispense mode: Type 2 Type 2

Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor	Default dilution
Standart	1.5				1:1.00	<input checked="" type="radio"/>
1:4	20.0	1.5	80		1:4.00	<input type="radio"/>
					2:1.00	<input type="radio"/>

Assay of

Done Cancel

?

# Имуноглобулин М, Ig-M

ImmunoTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: IgM Assay number: \*\*\* Date: 20.02.2020  
Calibration method: Spline Time: 17:14  
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator set: TruCalProt  
Replicates: 2 [Range 1 - 3]

Calibrator level: Blank: Water  
Cal 1: TruCalProt1  
Cal 2: TruCalProt2  
Cal 3: TruCalProt3  
Cal 4: TruCalProt4  
Cal 5: TruCalProt5  
Cal 6: None

Concentration: 0.0000

Assay of

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: IgM Assay number: \*\*\* Date: 20.02.2020  
Calibration method: Spline Time: 17:14  
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator: TruCalProt

Calibrator level	Sample	Diluted sample	Diluent	Water
Blank: Water	1.5			
Cal 1: TruCalProt1	1.5			
Cal 2: TruCalProt2	1.5			
Cal 3: TruCalProt3	1.5			
Cal 4: TruCalProt4	1.5			
Cal 5: TruCalProt 5	1.5			
Cal 6: None				

Assay of

Done Cancel

?

# Имуноглобулин М, Ig-M

ImmunoTD (Иммунотурбидиметрический)



## Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: IgM Assay number: \*\*\* Date: 20.02.2020  
Calibration method: Spline Time: 17:14  
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibration intervals:

Full interval: 0 (hours)  
Adjust interval: 0 (hours)

Calibration type:

Adjust type: 2 Point  
Adjust level: Cal 1  
Default ordering type: Full

Assay of

Done  
Cancel

?

Единицы измерения - г/л; знаков после запятой 2 (**Configure result units**). Диапазон метода: 0.25 - 8 г/л (**Configure assay parameters\ Results**).

Референсный интервал: 0.4-2.3 г/л (**Configure assay parameters\ Results**).  
При наличии ЛИС, не обязательно.

Диапазон нормальных значений указан в соответствии с рекомендациями производителя реагентов.

Калибратор TruCal Protein, набор из 5 уровней, нулевой - вода (или физраствор).

Контроль по TruLab Protein уровень 1 и уровень 2.