

Ферритин (DiaSys)

ИммуноTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Ferrit Type: Photometric Version: 1
Number: *** Assay availability: Enabled Date: 18.02.2020
Run controls for onboard reagents by: Lot Time: 11:01:08
Operator: ADMIN

Reaction definition Reagent / Sample Validity checks

Reaction mode: Rate up

Primary Secondary Read times
Wavelength: 604 / None Main: 19 -- 27
Last required read: 33 Flex: --
Absorbance range: -- Color correction: --
Sample blank type: None

Assay of

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Ferrit Type: Photometric Version: 1
Number: *** Assay availability: Enabled Date: 20.02.2020
Run controls for onboard reagents by: Lot Time: 17:14:07
Operator: FSE

Reaction definition Reagent / Sample Validity checks

Reagent: Ferrit
Diluent: Saline
Diluent dispense mode: Type 0

R1 R2
Reagent volume: 107 54
Water volume: --
Dispense mode: Type 2 Type 2

| Dilution name | Sample | Diluted sample | Diluent | Water | Dilution factor | Default dilution |
|---------------|--------|----------------|---------|-------|-----------------|----------------------------------|
| Standart | 5.4 | | | | = 1:1.00 | <input checked="" type="radio"/> |
| 1:5 | 20.0 | 5.4 | 80 | | = 1:5.00 | <input type="radio"/> |
| 1:2 | 2.7 | | | | = 1:2.00 | <input type="radio"/> |

Assay of

Done Cancel

?

Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Ferrit Assay number: *** Date: 20.02.2020
 Calibration method: Spline Time: 17:14
 Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator set: TruCalFrn
 Replicates: 2 [Range 1 - 3]

Calibrator level: Blank: Water
 Cal 1: TruCalFrn1
 Cal 2: TruCalFrn2
 Cal 3: TruCalFrn3
 Cal 4: TruCalFrn4
 Cal 5: None
 Cal 6: None

Concentration: 0.0000

Assay of

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Ferrit Assay number: *** Date: 20.02.2020
 Calibration method: Spline Time: 17:14
 Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator:

| Calibrator level | Sample | Diluted sample | Diluent | Water |
|-------------------|--------|----------------|---------|-------|
| Blank: Water | 5.4 | | | |
| Cal 1: TruCalFrn1 | 5.4 | | | |
| Cal 2: TruCalFrn2 | 5.4 | | | |
| Cal 3: TruCalFrn3 | 5.4 | | | |
| Cal 4: TruCalFrn4 | 5.4 | | | |
| Cal 5: None | | | | |
| Cal 6: None | | | | |

Assay of

Done Cancel

?

Ферритин (DiaSys)

ImmunoTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Ferrit Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibration intervals:

Full interval: 0 (hours)
Adjust interval: 0 (hours)

Calibration type:

Adjust type: 2 Point
Adjust level: Cal 1
Default ordering type: Full

Assay of

Done
Cancel

?

Единицы измерения - мкг/л; знаков после запятой 0 (Configure result units).

Диапазон метода 5-1000 мкг/л (Configure assay parameters\ Results).

Референсный интервал 15-400 мкг/л (Configure assay parameters\ Results).
При наличии ЛИС не обязательно.

Диапазон нормальных значений указан в соответствии с рекомендациями производителя реагентов.

Калибратор TruCal Ferritin, набор из 4-х уровней, нулевой - вода (или физраствор).

Контроль по TruLab Protein уровень 1 и уровень 2.