

Имуноглобулин G, Ig-G

ImmunoTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

☒ General ☐ Calibration ☐ SmartWash ☐ Results ☐ Interpretation

Assay: Type: Photometric Version: 1
Number: Assay availability: Date: 18.02.2020
Run controls for onboard reagents by: Time: 11:01:08
Operator: ADMIN

☒ Reaction definition ☐ Reagent / Sample ☐ Validity checks

Reaction mode:

Primary Secondary Read times
Wavelength: / Main: --
Last required read:
Absorbance range: -- Color correction: --
Sample blank type: Blank: --

Assay 13 of 54

Done Cancel

?

Configure assay parameters...

☒ General ☐ Calibration ☐ SmartWash ☐ Results ☐ Interpretation

Assay: Type: Photometric Version: 1
Number: Assay availability: Date: 20.02.2020
Run controls for onboard reagents by: Time: 17:14:07
Operator: FSE

☐ Reaction definition ☒ Reagent / Sample ☐ Validity checks

Reagent:
Diluent:
Diluent dispense mode:

R1 R2
Reagent volume:
Water volume:
Dispense mode:

Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor	Default dilution
Standart	1.5				1:1.00	<input checked="" type="radio"/>
1:4	20.0	1.5	80		1:4.00	<input type="radio"/>
					2:1.00	<input type="radio"/>

Assay of

Done Cancel

?

Имуноглобулин G, Ig-G

ИммуноTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: IgG Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator set: TruCalProt
Replicates: 2 [Range 1 - 3]

Calibrator level: Blank: Water
Cal 1: TruCalProt1
Cal 2: TruCalProt2
Cal 3: TruCalProt3
Cal 4: TruCalProt4
Cal 5: TruCalProt5
Cal 6: None

Concentration: 0.0000

Assay of

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: IgG Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator: TruCalProt

Calibrator level	Sample	Diluted sample	Diluent	Water
Blank: Water	1.5			
Cal 1: TruCalProt1	1.5			
Cal 2: TruCalProt2	1.5			
Cal 3: TruCalProt3	1.5			
Cal 4: TruCalProt4	1.5			
Cal 5: TruCalProt 5	1.5			
Cal 6: None				

Assay of

Done Cancel

?

Имуноглобулин G, Ig-G

ImmunoTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: IgG Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibration intervals:

Full interval: 0 (hours)
Adjust interval: 0 (hours)

Calibration type:

Adjust type: 2 Point
Adjust level: Cal 1
Default ordering type: Full

Assay of

Done
Cancel

?

Единицы измерения - г/л; знаков после запятой 2 (**Configure result units**).

Диапазон метода: 1.5 - 35 г/л (**Configure assay parameters\ Results**).

Рефереренсный интервал: 7-16 г/л (**Configure assay parameters\ Results**).

При наличии ЛИС, не обязательно.

Диапазон нормальных значений указан в соответствии с рекомендациями производителя реагентов.

Калибратор TruCal Protein, набор из 5 уровней, нулевой - вода (или физраствор).

Контроль по TruLab Protein уровень 1 и уровень 2.