




Certificate of Traceability

Manufacturer Certification

This certifies that **Shenzhen Lifotronic Technology Co., Ltd.**, using **H9 Hemoglobin Analyzer (HPLC)** has participated in and successfully completed the NGSP certification for manufacturers and is traceable to the **Diabetes Control and Complications Trial** Reference method. The comparison was performed with: **Zhongshan Hospital ASRL#3**

The system evaluated was:

Instrument: H9 Hemoglobin Analyzer (HPLC)	Calibrator Lot: 9SX3230020	Column Lot: SX23C43614
Reagent Lot: 9SX3230010	Calibrator Assigned Values: 5.60%, 10.30%	

Date of Certification: February 1, 2024

Certification Expires: February 1, 2025

NGSP Steering Committee Chair

NGSP Network Coordinator

SRL director/ supervisor